DD Scanbody 2.0

DE Gebrauchsanweisung, **EN** Instructions for use, **US** Instructions for use, **FR** Notice d'utilisation, **IT** Istruzioni per l'uso, **ES** Instrucciones de uso, **NL** Gebruiksaanwijzing, **EL** Οδηγίες χρήσης για το, **PT** Instruções de utilização, **TR** Kullanım talimatı, **RU** Инструкция по применению, **SV** Bruksanvisning, **ET** Kasutusjuhend, **DA** Brugsanvisning, **NO** Bruksanvisning, **FI** Käyttöohje, **IS** Notkunarleiðbeiningar, **PL** Instrukcja używania, **CS** Návod k použití skenovacího tělíska, **HU** Használati útmutató, **RO** Instrucțiuni de utilizare, **SK** Návod na použitie, **BG** Инструкция за употреба, **SR** Упутство за употребу, **FA** راهنمای استفاده از, **AR** دليل استخدام, **ZH** 使用说明书, **KO** 사용 설명서, **HR** Upute za uporabu, **GA** Treoracha úsáide, **LV** lietošanas pamācība, **LT** Naudojimo instrukcija, **SL** Navodila za uporabo, **MT** Struzzjonijiet għall-użu, **JA** 取扱説明書, **BS** Uputstvo za upotrebu

DE	DEUTSCH	Gebrauchsanweisung	4
			5
EN	ENGLISH	Instructions for use	6
			7
US	UNITED STATES (FDA)	Instructions for use	8
			9
FR	FRANÇAIS	Notice d'utilisation	10
			11
IT	ITALIANO	Istruzioni per l'uso	12
			13
ES	ESPAÑOL	Instrucciones de uso	14
			15
NL	NEDERLANDS	Gebruiksaanwijzing	16
			17
EL	ΕΛΛΗΝΙΚΑ	Οδηγίες χρήσης για το	18
			19
PT	PORTUGUÊS	Instruções de utilização	20
			21
TR	TÜRKÇE	Kullanım talimatı	22
			23
RU	РУССКИЙ	Инструкция по применению	24
			25
SV	SVENSK	Bruksanvisning	26
			27
ET	EESTI	Kasutusjuhend	28
			29
DA	DANSK	Brugsanvisning	30
			31
NO	NORSK	Bruksanvisning	32
			33
FI	SUOMI	Käyttöohje	34
			35
IS	ÍSLENSKA	Notkunarleiðbeiningar	36
			37
PL	POLSKI	Instrukcja używania	38
			39

CS	ČEŠTINA	Návod k použití skenovacího tělíska	40
			41
HU	MAGYA	Használati útmutató	42
			43
RO	ROMÂNĂ	Instrucțiuni de utilizare	44
			45
SK	SLOVENSKY	Návod na použitie	46
			47
BG	БЪЛГАРСКИ	Инструкция за употреба	48
			49
SR	СРПСКИ	Упутство за употребу	50
			51
FA	فارسی	راهنمای استفاده از	52
			53
AR	عربي	دليل استخدام	54
			55
ZH	中文	使用说明书	56
			57
KO	한국어	사용 설명서	58
			59
HR	HRVATSKI	Upute za uporabu	60
			61
GA	GAEILGE	Treoracha úsáide	62
			63
LV	LATVIEŠU VALODA	lietošanas pamācība	64
			65
LT	LIETUVIŲ KALBA	Naudojimo instrukcija	66
			67
SL	SLOVENŠČINA	Navodila za uporabo	68
			69
MT	MALTI	Struzzjonijiet għall-użu	70
			71
JA	日本語	取扱説明書	72
			73
BS	BOSANSKI	Uputstvo za upotrebu	74
			75

Inhalt:

- 1 x DD Scanbody 2.0
- 1 x DD Lab Screw, blau
- 1 x DD Implant Screw, silberfarben

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Contents:

- 1 x DD Scanbody 2.0
- 1 x DD Lab Screw, blue
- 1 x DD Implant Screw, silver

Product description:

The Dental Direkt Scanbody 2.0 is used for digital acquisition of the implant position relative to the residual dentition and soft tissue. The scan bodies are intended for optical acquisition and are screwed into the implant as a placeholder for intraoral scans or into the laboratory analog in the dental model. They are made of Ti6Al4V (medical grade 5) for single use and are supplied non-sterile, each with a separately packaged Implant and Lab Screw.

Caution: Dental Direkt scan bodies must only be used in combination with the DD Implant Library.

Intended purpose:

DD Solid Connect® implant components are used to support prosthetic restorations in combination with endosseous dental implants in the upper and/or lower jaw.

Indication:

For digital acquisition of the implant position on the model or for short-term use in the mouth, relative to the residual dentition and soft tissue.

Please note: DD Solid Connect® implant components from a particular series can only be combined with the corresponding compatible implant system. This also applies to the platform diameters.

Contraindication:

- Intolerance to the constituents
- Restorations with an angulation correction of more than 30° to the implant axis

Please note: Observe the contraindications specified by the original implant manufacturer used in each case.

Intended users:

DD Solid Connect® implant components should only be used by dentists and dental technicians in compliance with the procedures in the instructions for use.

Intended patient group:

Fixed restorations using Dental Direkt Solid Connect® implant components are suitable for the permanent dentition in adult patients of any gender and nationality.

Instructions for use in the laboratory:

The Scanbody 2.0 is supplied with a Lab Screw for use on the model and with a silver Implant Screw for use in the patient's mouth. Before scanning, the surface must be clean and free of grease.

The Scanbody 2.0 should be checked before scanning to make sure it does not show any signs of damage. Error-free scanning and exact reproduction of the implant position is only assured if the Scanbody is intact and used correctly.

Caution: Only use the Lab Screw in the model (max. torque: 10 Ncm).

Instructions for use in the dental practice:

The Scanbody 2.0 is inserted with the silver Implant Screw into the cleaned endosseous dental implant located in the patient's mouth.

Caution: The Lab Screw must not be used for placement of the implant components in the patient's mouth.

Caution: The Implant Screw must be tightened to the torque specified for the system (see table).

Material:

Scanbody 2.0, DD Implant Screw and DD Lab Screw are all made of Ti6Al4V (medical grade 5) ASTM 136 ELI.

Cleaning, disinfection and sterilization:

The Scanbody 2.0 and accompanying screws are cleaned at Dental Direkt, but not sterilized, before being packaged. Before they are used on the patient, the implant components must be cleaned, disinfected and sterilized following the procedure described by Dental Direkt.

Cleaning/disinfection:

Automatic cleaning in a suitable washer-disinfector with an appropriate cleaning program for medical instruments (e.g. "Vario TD Program", Miele).

OR

The following manual cleaning procedure:

- Rinse with cold tap water for 1 minute
- Rinse the screw channel using a disposable syringe (min. capacity 10 ml)
- Clean in the ultrasonic cleaner at 60°C (140°F) +/- 2°C, e.g. using Dr. Weigert neodisher® LM 2 for 10 minutes at a concentration of 5%
- Neutralize with cold distilled water for 2 minutes
- Rinse in the ultrasonic cleaner using distilled water at 60°C (140°F) for 10 minutes

Caution: The Scanbody 2.0, DD Implant Screw and DD Lab Screw are suitable for single use only. Any reuse may lead to hazards due to infection.

Sterilization:

Double-seal products in sterilization film conforming to EN ISO 11607-1. Steam sterilize in a suitable sterilizer with 3x fractionated prevacuum:

134°C (273°F) for 5 minutes, drying time: 20 minutes.

Then store in sterile packaging.

Possible side effects and interactions:

Allergic reactions or skin irritations are possible in very rare cases. Different alloys in the oral cavity can lead to galvanic reactions when in direct contact.

Disposal:

In compliance with local regulations.

Warranty:

Our products undergo strict quality control and conform to the latest technical standards. Please refer to our current warranty terms for further details.

Please note: Any serious incident that occurs in connection with the product should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is located

We are continuously developing and enhancing our products, and therefore reserve the right to make changes. The latest version of the instructions for use can also be found on our website at:

www.dentaldirekt.de/en

This version replaces all previous versions.

Explanation of symbols:

Manufacturer

Date of manufacture

Batch code

Catalog number

Keep dry

Non-sterile

Do not use if package is damaged

Torque

Consult electronic instructions for use
www.dentaldirekt.de/IFU

Caution: Under U.S. federal law, the device may only be sold by or on behalf of a dentist.

Do not re-use

Medical device

Unique device identification

Dental Direkt Series	Compatible with: manufacturer	System	Size (mm)	Torque (Ncm)
DD AAS-Serie	DENTSPLY Implants®	Astra Tech®	3.5 / 4.0	25
			4.5 / 5.0	25
DD BSK-Serie	Bredent Medical	Sky®	3.5 – 5.5	25
DD CCA-Serie	Altatec	Camlog®	3.3	20
			3.8	20
			4.3	20
			5.0	20
DD CCO-Serie	Altatec	Conelog®	3.3	20
			3.8	20
			4.3	20
			5.0	20
DD DEV-Serie	DENTSPLY Implants®	Astra EV	3.0	25
			3.6	25
			4.2	25
			4.8	25
			5.4	25
DD DFX-Serie	DENTSPLY Implants®	Frialit®/Xive®	3.4	25
			3.8	25
			4.5	25
			5.5	25
DD MAR-Serie	MegaGen	AnyRidge®	3.5 – 8.0	30
DD MMI-Serie	Medentis Medical®	ICX	3.75 – 4.8	30
DD MSE-Serie	MIS Implants®	Seven®	3.3 NP	30
			3.75 / 4.2 SP	30
			5.0 WP	30
DD NBR-Serie	Nobel Biocare™	Brånemark™	3.5 NP	35
			4.1 RP	35
DD NNA-Serie	Nobel Biocare™	Nobel Active®	3.5 NP	25
			4.3 / 5.0 RP	35
DD NRE-Serie	Nobel Biocare™	Replace Select®	3.5 NP	35
			4.3 RP	35
			5.0 WP	35
DD OTS-Serie	Osstem Implant	TS	3.5 M	20
			4.0 – 6.0 R	20
DD SBO-Serie	Straumann®	Bone Level®	3.3 NC	35
			4.1 / 4.8 RC	35
DD SSY-Serie	Straumann®	Tissue Level	4.8 RN	35
			6.5 WN	35
DD TSP-Serie	Thommen Medical	SPI®	4.0	25
			4.5	25
			5.0	25
			6.0	25
DD ZBC-Serie	Zimmer Biomet 3i	Certain®	3.4	30
			4.1 / 5.0	30
DD ZTA-Serie	Zimmer Dental	Tapered Screw-Vent®	3.5	30
			4.5	30
			5.7	30

Contents:

- 1 x DD Scanbody 2.0
- 1 x DD Lab Screw, blue
- 1 x DD Implant Screw, silver

Indications for Use:

DD Solid Connect® Implant Components are used to support prosthetic restorations in combination with endosseous dental implants in the upper and/or lower jaw.

All digitally designed custom abutments for use with DD Solid Connect® Implant components are to be sent to a FDA registered milling center for manufacture, validated by Dental Direkt GmbH.

Device description:

The Dental Direkt Scanbody 2.0 is used for digital acquisition of the implant position relative to the residual dentition and soft tissue. The scan bodies are intended for optical acquisition and are screwed into the implant as a placeholder for intraoral scans or into the laboratory analog in the dental model. They are made of Ti6Al4V (medical grade 5) for single use and are supplied non-sterile, each with a separately packaged Implant and Lab Screw.

Caution: Dental Direkt scan bodies must only be used in combination with the DD Implant Library.

Indication:

DD Scanbody 2.0: For digital acquisition of the implant position on the model or for short-term use in the mouth, relative to the residual dentition and soft tissue.

DD Implant Screw: The implant screw (abutment screw) is used to fix the abutment to the implant (two-part implant system).

Material:

The DD Scanbody 2.0, DD Lab Screw and DD Implant Screw are made of Ti6Al4V (medical grade 5) ASTM 136 ELI.

Contraindication:

- Intolerance to the components contained
- Restorations with an angulation correction of more than 30° to the implant axis

Note: Observe the contraindications specified by the original implant manufacturer used in each case.

Intended users:

The DD Solid Connect® Implant Components may only be processed by trained, dental technical / dental medical personnel in compliance with the specifications in the Instructions for Use.

Intended patient group:

Fixed dentures made of Dental Direkt Solid Connect® Implant Components are suitable for the permanent dentition of adult patients of any gender and nationality.

Possible side effects and interactions:

Allergic reactions or skin irritations are possible in very rare cases. Different alloys in the oral cavity can lead to galvanic reactions when in direct contact.

Use in MR (Magnetic Resonance) environment:

The DD Solid Connect® Implant Components have not been evaluated for safety and compatibility in the MR environment. They have not been tested for heating, migration, or image artifact in the MR environment. The safety of DD Solid Connect® Implant Components in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

Cleaning, disinfecting and sterilizing:

The DD Scanbody 2.0 and the associated implant screws are cleaned at Dental Direkt, but packed non-sterile. After processing and before use on the patient, the implant components must be cleaned, disinfected and sterilized following the procedure described by Dental Direkt.

Cleaning / disinfecting:

Automatic cleaning in a washer-disinfector suitable for this application with an appropriate cleaning program for medical instruments (e.g. "Vario TD Program", Miele).

OR

the following manual cleaning procedure:

- Rinse with cold tap water for 1 minute
- Rinse the screw channel using a disposable syringe (min. capacity 10 ml)
- Clean at 60°C +/- 2°C, e.g. using Dr. Weigert neodisher® LM 2 for 10 minutes at a concentration of 5% in the ultrasonic cleaner
- Neutralize with cold distilled water for 2 minutes
- Rinse in the ultrasonic cleaner using distilled water at 60°C for 10 minutes

Sterilizing:

Double-seal the product in FDA cleared sterilization film conforming to ISO 11607-1 and sterilize according to validated parameters listed below.

Steam sterilize in a FDA cleared sterilizer with 3x fractionated prevacuum: 132°C (270°F) for 4 minutes, drying time: 20 minutes.

Then store in sterile packaging.

Caution: DD Scanbody 2.0, DD Implant Screw and DD Lab Screw are suitable for single use only. Reuse may result in hazards due to infections.

Instructions for use in the laboratory:

The Scanbody 2.0 is supplied with a Lab Screw for use on the model and with a silver Implant Screw for use in the patient's mouth. Before scanning, the surface must be clean and free of grease.

The Scanbody 2.0 should be checked before scanning to make sure it does not show any signs of damage. Error-free scanning and exact reproduction of the implant position is only assured if the Scanbody is intact and used correctly.

Caution: Only use the Lab Screw in the model (max. torque: 10 Ncm).

Instructions for use in the dental practice:

The Scanbody 2.0 is inserted with the silver Implant Screw into the cleaned endosseous dental implant located in the patient's mouth.

Caution: The Lab Screw must not be used for placement of the implant components in the patient's mouth.

Caution: The Implant Screw must be tightened to the torque specified for the system (see table).

Note: DD Solid Connect® Implant Components of the respective series can only be combined with the corresponding compatible implant system. This also applies to the corresponding platform diameters.

Disposal:

Dispose of waste according to applicable legislation.

Warranty:


Our products are subject to strict quality control and correspond to the current technical standards. Please refer to our current warranty terms and conditions for further details.

Note: Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

Our products are subject to continual development, so we reserve the right to make changes. You will find the latest version of the instructions for use on our homepage at: www.dentaldirekt.de/en.

This version replaces all previous.

Explanation of symbols:

-  Manufacturer
-  Date of manufacture
-  Batch code
-  Catalog number
-  Keep dry
-  Non-sterile
-  Do not use if package is damaged
-  Torque
-  Consult electronic instructions for use www.dentaldirekt.de/IFU
-  Caution: U.S. Federal law restricts this device to sale by or on the order of a dentist
-  Do not re-use
-  Medical Device
-  Unique Device Identification

Dental Direkt Series	Compatible with: manufacturer	System	Size (mm)	Torque (Ncm)
DD AAS-Serie	DENTSPLY Implants®	Astra Tech®	3.5 / 4.0	25
			4.5 / 5.0	25
DD CCA-Serie	Altatec	Camlog®	3.3	20
			3.8	20
			4.3	20
			5.0	20
DD DEV-Serie	DENTSPLY Implants®	Astra EV	3.0	25
			3.6	25
			4.2	25
			4.8	25
			5.4	25
DD DFX-Serie	DENTSPLY Implants®	Frialit®/Xive®	3.4	25
			3.8	25
			4.5	25
			5.5	25
DD NNA-Serie	Nobel Biocare™	Nobel Active®	3.5 NP	25
			4.3 / 5.0 RP	35
DD NRE-Serie	Nobel Biocare™	Replace Select®	3.5 NP	35
			4.3 RP	35
			5.0 WP	35
DD SBO-Serie	Straumann®	Bone Level®	3.3 NC	35
			4.1 / 4.8 RC	35
DD SSY-Serie	Straumann®	Tissue Level	4.8 RN	35
			6.5 WN	35
DD ZBC-Serie	Zimmer Biomet 3i	Certain®	3.4	30
			4.1 / 5.0	30
DD ZTA-Serie	Zimmer Dental	Tapered Screw-Vent®	3.5	30
			4.5	30
			5.7	30

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