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Symbolerklärungen

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Instruction Manual

Explanation of symbols

 Keep dry	 Use-by date	 Caution: Health hazard
 Consult instructions for use	 Temperature limit 18°C / 64°F 28°C / 82°F	 Caution: Serious health hazard
 Batch code	 Catalogue number	 Caution: Hazardous to the environment
 Keep away from sunlight		

1. Indications

DD exaprint is a light-curing dental resin for the production of dental working and demonstration models in 3D printers (DLP; SLA; LCD).

2. Contraindications

The material is not a medical device and should only be used for the intended indication and purpose of 3D printing.

3. Safety information

Please take note of the information in the latest version of the safety data sheet. Avoid direct skin contact with non-polymerized (liquid) material during processing. Wear protective gloves, protective goggles and a face mask to avoid irritation. If the product is not used in line with the instructions for use, the quality of the product will be impaired.

4. Handling and storage

Only use clean material from undamaged packaging. Store the material in tightly sealed containers out of sunlight at 18°C-28°C. **DD exaprint** should only be exposed to daylight for short periods during use otherwise premature polymerization of the material will occur. To avoid contamination between different printing materials, we recommend that you use separate material trays.

5. Guidelines for use

Please adhere to the following protocol in order to achieve optimal results:

Processing:

Make sure in advance that the 3D printer is in proper condition, clean and filled with sufficient material.

- Gently tilt the material container back and forth several times before use - do not shake to avoid bubbles from forming!
- Carefully fill the resin into the material tray and let it rest until no more bubbles are visible.
- Maintain a processing temperature of 23-35°C (depending on the ambient temperature).

Design conditions to be considered:

- Recommended layer thickness of the DLP print: 50 µm
- 3 mm wall thickness recommended for cavity position of the model.

After the printing process:

- We recommend a building platform drip off time of approx. 10 minutes after the printing process has been completed.
- Immediate post-processing of the print object necessary to avoid deformations. If immediate post-processing is not possible, leave the print objects on the building platform in the printer.
- Remove residual material from the material tray after printing - within 12h at the latest back into the original container or, if necessary, into a separate light-impermeable container.

Cleaning:

- Use an unheated ultrasonic bath filled with isopropanol (≥97%) for the cleaning process. Follow the operating instructions/instructions for use issued by the device manufacturer.

The print objects are cleaned in two steps:

Step 1 - Pre-cleaning process:

- **Optional:** Before the actual pre-cleaning process is carried out, the printing objects can be pre-treated by repeatedly immersing them in an isopropanol bath for 15 seconds. The print objects must be held on a support using forceps.
- 3 minutes in an ultrasonic isopropanol bath. The openings of the print object should face downwards. The bath for pre-cleaning can be used several times.
- However, the cleaning fluid should be changed whenever the objects do not come out as clean.

Step 2 - Main cleaning process:

- Clean print objects for a further 3 minutes in a fresh isopropanol ultrasonic bath.
- Do not use cleaning fluid more than once in this step in order to achieve a consistent quality.

The print objects should then be carefully dried with compressed air. If there are still resin residues on the print objects after drying, or if they appear from undercuts during drying, then please repeat the main cleaning process.

Post-exposure:

Post-polymerization has a significant influence on the surface quality and hardness of the printed objects. It is recommended that post-exposure take place 20 minutes after the main cleaning process. Ensure that print objects do not touch each other or overlap. Shadow formation impairs post-polymerization.

Optionally, support structures that get in the way can be carefully removed before post-exposure. The print object must then be cleaned again using isopropanol. However, in general it is not necessary to remove this.

In principle, different devices can be used for post-exposure. For optimal results with **DD exaprint** resins we recommend:

- **Otoflash G171** – Polymerization device with xenon flash light (Manufacturer: NK-Optik GmbH)

Parameters for Otoflash G171

Item no.	Product	Parameters	Protective atmosphere
D120	DD exaprint	1. 2000 flashes 2. Turn print object – 2 min. cooling stage with open lid 3. 2000 flashes	Not necessary – however, use of nitrogen N2 (purity grade 2.6) is recommended

Final processing:

DD exaprint can be finished and polished in the usual way. To achieve optimum results, ensure that the handpiece has a low contact pressure and a low speed.

Processing residues can be removed with compressed air. The final cleaning is carried out by evaporation. Excessive heat should be avoided when doing this.

6. Disposal

Dispose of liquid, dissolved (in isopropanol) or partially hardened resin residues in accordance with applicable official regulations. Must not go down the drain or be disposed of together with household waste.

*Our instructions for use and notes have been prepared for you to the best of our knowledge and are forwarded to you without obligation. It is your sole responsibility to test delivered products and their application. The quality of our products is always guaranteed according to our standards. We are continuously developing our products, and therefore reserve the right to make changes. The latest version of the instructions for use can be found on our website at: www.dentaldirekt.de/en **This version replaces all previous versions.***